

Finally.

Aesthetics that your patients will love.



Which would you choose?

Zirlux® is a strong, highly aesthetic alternative to metal that provides patients with the aesthetics and functionality they demand.



Patients choose Zirlux®!



Metal-Occlusal Crown



Gold Crown



Prescribe Zirlux® Today!

Zirlux® is indicated for full contour crowns, bridges, inlays, onlays, and traditional layered ceramic restorations.



For More Information or to Prescribe Zirlux®
Contact Your Zirlux® Certified Partner Today!




1-866-248-9657 | www.zirlux.com


The Zirlux® Universal All-Ceramic System enables you to provide patients with exceptional restorations that are aesthetic, strong, metal-free, and available for multiple restorative options. Zirlux® can be prescribed for translucent full contour monolithic zirconia, traditional layered zirconia, or pressed to zirconia restorations.

Zirlux® Benefits

- High translucency for increased aesthetics
- Suitable for supra gingival and feather edge margins
- Pre-shaded material produces consistent and predictable restorations
- Aesthetic alternative to metal
- Low wear to opposing enamel



Dr. Burgess, et al.
28 Wear of Ceramics and Enamel in artificial chewing simulator
AADR, Tampa, FL 2012



* Ardent, Inc R&D 2011

Preparation Guidelines

General preparation guidelines for all ceramic restorations should be followed: (Fig 1. and Fig. 2)


Preparation should follow the anatomy of the tooth, providing at least the minimum thickness required for the respective restoration. A 1.0 mm axial and occlusal reduction is considered ideal for full contour zirconia restorations. (Table 1)

A definitive finish line (shoulder with rounded internal line angles or chamfer) is recommended. Feather edge preparations are acceptable. (Fig. 3)

All sharp edges and line angles should be rounded.

| TABLE 1 Minimum Zirconia Thickness - Crowns | | | |
|--|--------|--|--------|
| Anterior | | Crowns | |
| Axial | 0.5 mm | | |
| Incisal/Occlusal | 0.5 mm | | |
| Connector Dimensions | n/a | | |
| Posterior | | Crowns | |
| Axial | 0.5 mm | | |
| Incisal/Occlusal | 0.5 mm | | |
| Connector Dimensions | n/a | | |
| Minimum Zirconia Thickness - Bridges | | | |
| | | ≥4 units with mix 2 connected pontics | |
| | | Cantilever with 1 pontic | |
| Axial | 0.5 mm | 0.7 mm | 0.7 mm |
| Incisal/Occlusal | 0.7 mm | 1.0 mm | 1.0 mm |
| Connector Dimensions | 7 mm² | 9 mm² | 12 mm² |
| Posterior | | | |
| Axial | 0.5 mm | 0.7 mm | 0.7 mm |
| Incisal/Occlusal | 0.7 mm | 1.0 mm | 1.0 mm |
| Connector Dimensions | 9 mm² | 12 mm² | 12 mm² |


Ideal Chamfer Margins



Feather Edge Margin



FIGURE 2



Cementation and Bonding Reference

Zirlux® restorations may be cemented or bonded using conventional cements and luting materials. Resin modified glass ionomers and resin cements are recommended over Zinc Oxide phosphate cements to optimize aesthetic results. Inlays, onlays, and veneers must be adhesively cemented (bonded).

Ceramir® Crown & Bridge **DOXA**

RelyX™ Luting Cement **3M ESPE**

GC Fuji Plus® **'GC'**

Panavia™ F2.0 **kuraray**

Multilink® Automix **ivoclar vivadent**
passion vision innovation

Refer to the respective manufacturer's instructions for use for proper application of any cementation/bonding materials.

Contact your Zirlux® certified laboratory or Henry Schein Dental to discuss available cements and luting materials.

